**IFU-OTH-0003**

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Revisão / Revision: 0  
Data / Date: 07/09/2023

### Instructions for Use

**Technical Name:** Orthodontic Bracket

**Product Name:** Orthometric Metal Brackets

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- Telemóvel: +351 967 417 654

### Product description

Self-ligating brackets are made of stainless steel, with a self-ligating system in a resilient material and an anatomical design that establishes a good fit with the dental elements, providing the patient with comfort, quality, safety, and reliability.

When fixed to the buccal surface of the tooth, the brackets aim to house the intraoral orthodontic arches and interact with other orthodontic accessories, helping to establish the mechanical-biological system, thus enabling the application of force to the teeth, which results in induced tooth movement, correcting malocclusion. In some cases, brackets can serve as supports for the use of cantilevers, which increase their versatility in tooth movement.

U-Clip self-ligating brackets are found in the interactive clip self-ligating system, in which interaction occurs from the .016"X .022" rectangular arch. For the Ultra-P self-ligating bracket system, a passive port self-ligating system is used, which does not actively interact with the archwires inserted into the groove.

### Intended use

Brackets are used when orthodontic mechanics are required to correct the various malocclusions, sagittal, vertical and transverse, inter-arch and intra-arch, taking into account the appropriate diagnosis in all cases. The following malocclusions can be treated using brackets:

#### I- Sagittal problems:

- 1) Angle Class I malocclusion;
- 2) Angle Class II malocclusion;
- 3) Angle Class III malocclusion;
- 4) Mild mandibular prognathism (in compensatory treatment) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide adequate occlusion;

5) Mild mandibular retrognathism (in compensatory treatment in adults, and in patients in the active growth phase, making up the mandibular propulsion system) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;

6) Mild maxillary protrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adequacy;

7) Mild maxillary retrusion (compensatory treatment in adults, and in growing patients after the orthopedic phase) and in severe cases can be used in pre-surgical orthodontic preparation and in the post-surgical phase to provide occlusion adjustment;

8) Anterior crossbite (by correcting the malocclusion through buccal inclination of the upper incisors and/or lingual inclination of the lower incisors);

II- Transverse problems:

1) Posterior crossbite (in the post-orthopedic expansion phase or as a compensatory treatment for the buccal inclination of postero-superior teeth) and also after orthognathic surgery has been performed and the post-surgical containment phase has elapsed.

III- Vertical problems:

1) Deep overbite;

2) Anterior open bite (in orthodontic or orthodontic-surgical treatment) and in severe cases they can be used in pre-surgical orthodontic preparation and in the post-surgical phase to adjust the occlusion.

IV- Inter-arch and intra-arch problems:

1) Dental crowding;

2) Midline deviations;

3) Dental positioning asymmetries.

#### **How to use:**

For proper bonding of brackets, the following sequence of procedures must be adopted and are summarized as:

1. prophylaxis;

2. acid etching of the enamel;

3. application of the primer followed by bonding the bracket to the tooth.

The direct and indirect bonding techniques are described below:

#### **I- Direct bonding**

1) Isolation of the working area with plastic cheek retractors and relative isolation with cotton rollers, keeping the high-powered saliva sucker on to avoid contamination of the bonding areas with saliva;

2) Delimiting the specific areas for bonding accessories;

3) Performing prophylaxis with pumice stone, water and Robson's brush to remove the enamel surface, the inducts that are attached, Nasmyth's membrane, the acquired film and other contaminating agents;

4) Wash the sanitized areas with plenty of water;

5) Condition the enamel surface, only in the area that will receive the bracket, using 37% phosphoric acid gel applied for 15 to 30 seconds. In deciduous teeth, conditioning should be for 30 to 60 seconds. This acid etching procedure aims to create micro porosities for the mechanical embedding of the bonding resin;

6) Thoroughly rinse the etched areas, keeping the high-powered suction device in position to avoid contamination of the enamel surfaces. If saliva contamination occurs, the etching procedure is repeated;

7) Drying the conditioned enamel surfaces with a jet of air free from contamination with water or oil;

8) Application of fluid resin or primers to the enamel surface of the bonding area with a brush or microbrush, without exerting pressure;

9) Light-curing with the light-curing device, at a time and ultraviolet emission defined by the professional, according to the treatment;

10) Keep the bracket fixed between the ends of the bonding clamp code 75000321, avoiding handling the accessories directly with the fingers, which could contaminate them;

11) Applying a thin layer of primer all over the base of the bracket, covering the entire retention area;

- 12) Positioning the bracket on the enamel surface using the bracket bonding forceps, applying light pressure;
- 13) Checking the occlusal distance or X;
- 14) Removal of excess resin; and
- 15) Light curing with the LEDX-T 2400 light curing device, in the fast Ortho function, for 3 seconds per bracket.

#### II- Indirect bonding

- 1) Selection of the brackets of interest to the professional;
- 2) Choosing the trays;
- 3) Prophylaxis of the upper and lower dental arches to eliminate the presence of any contaminating physical agent from the enamel surfaces;
- 4) Molding the upper and lower dental arches with irreversible hydrocolloid (alginate);
- 5) Obtaining type IV stone plaster models or digital models;
- 6) Delimitation of the long axes of the dental elements buccally and lingually, and demarcation of the EV points and X bonding heights on the plaster models;
- 7) Applying plaster insulation to the plaster surfaces in the places that will receive the resin brackets;
- 8) Leave to dry for 30 seconds;
- 9) Bonding the attachments with gold quality orthodontic resin, removing excess resin and checking the positioning of the attachments;
- 10) Light-curing with the LEDX-T 2400 light-curing appliance, in the fast Ortho function, for 3 seconds per bracket; and
- 11) Obtaining the transfer trays.

Direct and indirect bonding techniques may vary depending on the professional who performs them. Due to the large number of steps involved in the bracket bonding process, it is very difficult to control the moisture absorbed in the process, which can lead to contamination of the enamel surface and, consequently, failure to adhere to the tooth. The presence of moisture and/or contamination during the bonding process is the biggest cause of bracket failure, delaying and increasing treatment costs.

Brackets can undergo changes to their surface structure when exposed to the adverse conditions present in the oral environment. In the oral cavity, metal biodegradation usually occurs through electrochemical disruption. Although electrolytic attack seems to be the major cause of disruption, bacteria and their cellular degradation products, and selective interactions with gases such as oxygen and carbon dioxide can contribute to the degradation of materials in the oral cavity. Organic acids and enzymes in particular can also affect the various metals and the pH of the environment where the devices are installed, exerting a major influence on the corrosive process.

The medical product in question may only be fitted to a patient by a dental surgeon specializing in orthodontics.

#### **Bracket removal**

Once the treatment objectives have been achieved and the finishing touches have been made, the brackets are removed with the adhesives without damaging the tooth enamel structure. It is recommended that suitable, high-quality instruments are used to avoid discomfort for the patient and the risk of injury to soft tissue and dental structures.

The technical removal procedures involve the following steps:

- 1) Removal of excess resin that may be surrounding the base of the brackets, however, if the bonding has been carried out with technical care, the minimum possible amount of resinous material around the orthodontic accessory is observed;
- 2) Use of How straight pliers (number 110) or Weingart pliers (number 120), both of which are used to perform mesiodistal compression of the bracket fins, which causes deformation of the accessory and its release from the enamel surface by breaking the bracket-resin interface. The instruments used at this stage will depend on the professional's familiarity, with some using 346 pliers or metal or plastic "pistol grip" pliers.
- 3) Remove the adhesive using scraper pliers with Widia or titanium tips. Excess resin can also be removed with a pair of pliers;

4) Removal of composite residues with a multilaminar conical trunk drill made of tungsten carbide, adapted to a 9-blade low rotation contra-angle, CB 27, code 75.80.1004. You can also use the high-rotation 24-blade drill, code CF 375R, code 75.80.1005. These drills have inverted blades to preserve tooth enamel; and  
5) After removing the resin residue, the surface of the enamel should be finished with a Shofu Ceramisté polishing rubber. The final polishing is carried out with rubber cups and Spain white in water or cloth wheels in a contra-angle and diamond paste.

**Composition:**

U-Clip and Ultra-P Bracket Body - 17-4 PH Steel (Designation UNS S17400).

**STORAGE CONDITIONS:**

In intact packaging and, where applicable, after opening.

This product does not require temperature, humidity or light control during storage. The recommendation is to keep it in a clean, dry place.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

**WARNINGS AND SYMBOLOGY:**

It is essential that only the specialist installs the product on the patient, since in addition to the clinical examination, the professional must have the necessary knowledge to make a systematic description of the patient's dentofacial features, list the priority problems, describe possible solutions for each of the problems, considering the possible favorable or unfavorable interactions for each of the proposed attempts.

It is also the orthodontist's duty to determine the maximum expected gain from orthodontic therapy and the minimum loss, or to calculate the benefit versus risk ratio for each patient, with both risk and benefit including psychosomatic factors in order to determine the best management for each patient.

Patients should be instructed in oral hygiene, maintenance and food care. The following should be avoided: ice, coconut, hard and/or sticky sweets, caramels, nougats, jelly beans, green corn, popcorn and sugar cane. Biting hard objects such as pens and pencils has the same effect. Chewing gum and toffee often loosen the bands or rings that have been cemented on.

The high concentration of bacteria and fungi in the mouth can accelerate the corrosion of metal orthodontic appliances. The pH of the oral environment also has a major influence on the corrosive process, so a diet rich in acidic foods, such as excessive intake of soft drinks with acidic pHs, should be avoided as much as possible during treatment.

**Contraindications**


Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment conducive to the demineralization of tooth enamel.

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Revisão / Revision: 0  
Data / Date: 07/09/2023

	Do not re-use		Authorized representative in the European Community
	Consult instructions for use		European conformity
	Medical device		Caution
	Catalogue number		Non-sterile
	Batch code		Contains Nickel
	Manufacturer		Date of manufacture
	Use-by date		Keep away from sunlight
	Keep dry		Do not use if package is damaged